
Software Release Note

IMS 6f-E7, IMS 7f, IMS 7f-Auto, IMS 7f-GEO

Release 3.6 Patch #8

Release Date

Date	Release	Applicability
September 3, 2025	3.6.8	This software release applies to <ul style="list-style-type: none">Correction R1854, 1856, 1871

Modifications in Software Release 3.6.8 for IMS 6f-E7, IMS 7f, IMS 7f-Auto, IMS 7f-GEO

1. List of Bugs Fixed

ID	Application	Description
1854	Camera	Fixes a bug in 'field of view calibration' associated with Allied Vision Alvium G1-507C cameras.
1856	Actinis	Fixes a bug in 'Intro application' relating to the display of 'intro camera'.
1871	Acquisition	Fixes a bug with lost CKB data. Previously when restarting a 'CheckerBoard' acquisition and clicking 'Cancel' in the 'Save As window', the data from the previous acquisition was lost. This has been corrected.